CE 0086 Syringe
CE 0197 Needle

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GB Instruction leaflet

Description
Teosyal® Kiss is a transparent, non pyrogenic, sterile, cross-linked, hyaluronic acid, visco-elastic gel of non-animal origin. Each box contains one or two pre-filled syringes of Teosyal® Kiss. For each syringe, the box contains two sterile 27G1/2" needles and two traceability labels (one to be given to the patient, and one to be kept by the doctor in the patient's file). The volume of each syringe is shown on the packaging carton as well as on each syringe.

Composition:
Cross-linked hyaluronic acid25 mg
Phosphate buffer pH 7.3qs ad 1 ml

Indications
The therapeutic purposes of Teosyal® products are:
- to modify the anatomy of aged skin: restoration of volumes, filling of skin depression and wrinkle, restoration of the skin hydration, regulation of the Feuchthaltigkeits der Haut
- reconstructive surgery: filling of depressions due to scars, reconstruction of volumes lost by lipotatrafy.
Teosyal® Kiss is indicated for increasing lip volume and for the correction of moderate to deep lines at the corners of the mouth.

Mode of action
Teosyal® Kiss is injected into the mucosa for the treatment of lips and into the deep dermis for the correction of cutaneous depressions. Teosyal® Kiss thus creates a volume that corrects wrinkles and restores the lip volume. Teosyal® Kiss is biodegradable and is absorbed slowly over time. According to the depth of the wrinkle or the volume of the lips to be treated, normally one to three treatment sessions are necessary to obtain optimal level of correction. Periodic touch-up sessions following the treatment maintain the desired level of correction.

Contra-indications
Teosyal® Kiss must not be used:
- for the correction of superficial wrinkles.
- for injections other than intradermal
- for non cutaneous or intradermal injections
- in combination with peeling, laser treatment or ultrasound-based
- if the patient has cutaneous disorders, inflammation or infection at the treatment site or near to this site.
- in the case of patients having a known hypersensitivity to hyaluronic acid, with a history of severe allergy or anaphylactic shock,
- in case of patients with autoimmune diseases,
- Due to possible interactions with other filling implants, which have not been researched, it is inadvisable to inject Teosyal® Kiss into sites in the presence of other filling implants.
- in pregnancy, breast-feeding mothers, or in children
Do not inject into blood vessels

Posology and method of administration
Teosyal® Kiss is only to be carried out by practitioners trained in injection techniques for products designed to correct wrinkles and increase lip volume. Before beginning treatment, patients must be questioned regarding their past medical history and informed regarding the foreseeable outcome of treatment and of potential undesirable effects. Correction sites must be properly disinfected. It is advised to perform injections using the needles provided in the box. Teosyal® Kiss is injected slowly into the mid to deep dermis to correct wrinkles or into the mucosa to treat lips, using a linear threading injection technique. If the injection is made too deeply, i.e., into subcutaneous tissue, the correction will not have the sought after duration of action. An injection which is too deep, will encounter subcutaneous tissue which will not give rise to resistance during injection of the product, in contrast to the dermis. If the colour of the needle is visible under the skin during injection, this indicates that the injection will be too superficial, which it is also advisable to avoid, as the result of the correction might be irregular. The injection volume is dependent upon the correction required. It is however recommended not to exceed 2 ml per session. It is important not to overcorrect. Massage the treated sites carefully in order to ensure a uniform distribution of the product at the corrected sites.

Side effects
The practitioner must inform the patient that there are potential side effects related to the implantation of this device, occurring immediately following injection or after a delay. These include, among others:
- Local manifestations: inflammatory reactions (erythema, oedema, pain at the point of injection), haematomas, itching, temporary loss of sensitivity around the injected area, dyschromia, abscesses, indurations, nodules (possibly granulomas), skin necrosis, migration of the implant.
- General manifestations: immediate hypersensitivity up to anaphylactic shock.
The patient should inform the practitioner of any secondary effect other than those described above or any effect persisting beyond one week.
The practitioner will, in turn, inform the product retailer within the shortest time possible.

Assembly the needle on the syringe
For optimum use of Teosyal® products, it is important that the needle be assembly on the syringe following the 4 steps detailed in the diagrams (see figures 1 to 5). Stop the injection and change the needle if you feel an obstruction or pressure during the injection.

Precautions for use
Verify the expiry date and the integrity of the packaging before use. Do not use if the expiry date is exceeded or if the packaging is damaged. You should warn the patient against taking high-dose Vitamin E, aspirin, anti-inflammatories or anti-coagulants, during the week before the injections session. Advise the patient not to use make-up during the 12 hours following injection and to avoid extreme temperatures (intense cold, sauna,...) during the week following treatment. At the end of the treatment session, it is essential to discard all remaining unused product. This device is intended for single use only. Sterility is not guaranteed in the event of re-use. The gel can deteriorate, reducing its efficiency and the gel can dry out, hindering its extrusion through a needle. Dispose of used needles in appropriate containers.

Storage conditions
Store between 2°C and 25°C, away from direct sunlight.
Make sure there are no visible signs of damage to the packaging before use.

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